




Emergo Europe  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands



Syvek is a registered trademark. US Patent Nos.  
5,622,834, 5,623,064, 5,624,679, 5,635,493,  
5,686,115, 5,846,952, 5,858,350, 6,063,911,  
6,599,720, 6,610,668, 6,630,459, 6,649,599,  
6,686,342, 6,743,783, 6,864,245, 7,041,657,  
7,115,588, 7,285,266, 7,931,637, 8,152,750,  
8,481,512, 8,486,033

MPS 500-62 REV-03 10/2015

**Syvek<sup>®</sup>Excel**  
PATCH



400-16



**MPT** MARINE POLYMER  
TECHNOLOGIES

Innovative Medical Solutions

107 Water St. Danvers, MA 01923  
1-800-666-2560

Made in USA

## DESCRIPTION

Syvek®ExcelPatch is a soft, white, sterile, lyophilized pad of a cellulosic polymer (poly-N-acetylglucosamine) isolated from microalgae, and attached to medical grade foam backing.

## INDICATIONS

Syvek®ExcelPatch is intended for use under the direction of a healthcare professional.

Syvek®ExcelPatch promotes the rapid control of bleeding in patients following hemodialysis and in patients on anticoagulation therapy.

Syvek®ExcelPatch is indicated for use in the local management of bleeding wounds such as vascular access site, percutaneous catheters or tubes and surgical debridement.

## PRECAUTIONS

There are no known contraindications to the use of Syvek®ExcelPatch. Longer compression may be necessary for hypertensive or obese patients.

### **Not Made with Natural Rubber Latex**

Store at controlled room temperature  
15°C to 30°C (59°F to 86°F).

Sterilized with radiation.

Do not resterilize and/or reuse.

One-time use only.

## Syvek®ExcelPatch APPLICATION

1. Open Syvek®ExcelPatch tray into sterile field.
2. Hold occlusive pressure proximal to the puncture site and remove the sheath.
3. Allow 1-2 drops of blood from site.
4. Immediately place Syvek®ExcelPatch (beige side AWAY from skin) over the puncture site.
5. Apply firm pressure over the puncture site.
6. Slowly decrease amount of proximal pressure after 4 to 6 minutes (3 - 4 minutes for diagnostic procedures) while continuing to apply firm compression over the puncture site.
7. Hold firm pressure on the puncture site until hemostasis is achieved, then slowly release.
8. Assess the site without disturbing the Syvek®ExcelPatch, then cover with a transparent adhesive dressing.

## REMOVE Syvek®Excel Patch after 24 Hours:

1. Remove the transparent adhesive dressing.
2. Soak Syvek®ExcelPatch with water.
3. Gently peel off Syvek®ExcelPatch.

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner.