Latex-Free, Powder-Free Examination Gloves
Gants d'examen sans latex et sans poudre

Ansell

www.ansellcanada.ca

Selecting Latex-Free Exam Gloves

Choosing the right glove is one of the most important decisions for healthcare workers today. The two primary considerations should always be barrier protection and allergen content. Because of allergy concerns, a growing number of healthcare facilities have eliminated or limited latex products, including gloves. The alternatives for latex examination gloves are latex-free materials including nitrile, neoprene, and vinyl.

When choosing an exam glove, the level of exposure to risk the patient-care activities require is a major selection criteria. Procedures that involve exposure to infectious material require a glove that provides appropriate barrier protection. It is important the gloves also fit well around the wrist. The Centers for Disease Control and Prevention (CDC) recommends that gloves fit the user's hand comfortably and should not be too loose or too tight. The following is a summary of the strengths and weaknesses of each latex-free gloving material.

Nitrile

Nitrile is extremely strong and offers puncture and chemical resistance superior to all glove films. Its tensile strength is also typically better than that of other gloving materials and its elasticity is extremely high, allowing the gloves to conform to the shape of the wearer's hands, providing good comfort and fit. There are no natural rubber latex (NRL) proteins in nitrile, so there is no chance of Type I latex allergy with use. Nitrile gloves are recommended as an alternative to NRL products for conditions where high strength and chemical protection are required, and for individuals with NRL allergies.

Neoprene (Polychloroprene)

Neoprene offers a similar fit, feel and barrier protection to latex. It's strong and provides resistance to many chemicals (acids, alcohols, caustics, detergents, and ketones). It is comfortable to wear for long periods, and its elasticity is close to that of latex with very high memory. The film is able to retain its original shape and is somewhat puncture resistant. Neoprene gloves are recommended for individuals with latex allergies.

Polyvinyl Chloride (PVC)

Vinyl works well for selected tasks involving limited patient contact. It may not provide a snug fit on the hand, especially around the wrist. Vinyl should always be chosen with caution because it provides reduced barrier protection when compared to other non-latex materials such as nitrile and neoprene. Vinyl offers poor resistance to chemicals such as the alcohols used in swabbing, and should not be used when handling chemotherapy drugs. In addition, its use should be limited to short procedures in low-risk situations that do not involve exposure to blood, body fluids, or chemicals. While thought to be less expensive, a Tufts University study concluded that the economic advantages of vinyl are overstated and that substituting vinyl with safer alternatives is cost effective and practical.

Safety and Protection

Infection Control

Viral Penetration. Although not required by the FDA, Ansell's premium line of gloves is tested to ensure protection against viral penetration. Our gloves meet the viral penetration requirements according to ASTM F1671.

Cuff. Beaded.

Skin Protection

Powder-Free. Helps eliminate powder-related complications and sensitivity.

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Latex-Free, Powder-Free Examination Gloves

In addition to the overall benefits that Ansell latex-free exam gloves offer, each glove provides unique features.

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Proud Partner

Susan G. Komen for the Cure®

Ansell is proud to partner with Susan G. Komen for the Cure® to unite in the fight against breast cancer.

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MICRO-TOUCH® NITRAFREE™

- Nitrile
- Chemical accelerator-free*
- Reduces the risk of Type IV allergy (chemical)
- Pink color
- Textured fingertips
- Tested for use with chemotherapy drugs ASTM D6978
- Length: 240 mm (9.4")

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* Does not contain sulphur-based accelerators such as thiurams, carbamates, thiazoles, guanidine or thioureas.

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MICRO-TOUCH® NITRA-TOUCH®

- Nitrile
- Blue color
- Textured fingertips
- Excellent chemical and puncture resistance
- Tested for use with chemotherapy drugs ASTM D6978
- Length: 245 mm (9.6")
- CGSB certified

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MICRO-TOUCH® NITRILE

- Nitrile
- Blue color
- Textured finish
- Excellent chemical and puncture resistance
- Tested for use with chemotherapy drugs ASTM D6978
- Length: 239 mm (9.4")

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MICRO-TOUCH® NITRILE E.P.®

- Nitrile
- Blue color
- Textured fingertips
- Excellent chemical and puncture resistance
- Tested for use with chemotherapy drugs ASTM D6978
- Length: 295 mm (11.6")

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NEOPRENE

MICRO-TOUCH® AFFINITY®

- Neoprene
- Green color
- Textured fingertips
- Easy donning due to polyurethane coating
- Tested for use with chemotherapy drugs
- Length: 239 mm (9.4")

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
VINYL

MICRO-TOUCH® ELITE®

- Stretch polyvinyl chloride
- Cream color
- Smooth finish
- Stretch for better fit and feel
- For low-risk applications
- Length: 241 mm (9.5")
- CGSB Certified

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
SYNSATION®

- Polyvinyl chloride
- Clear
- Smooth finish
- For low-risk applications
- Length: 241 mm (9.5")

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Barrier Protection and Risk Factors

There are significant differences in the durability of gloving materials, and it is necessary for healthcare workers to select gloves with the barrier durability appropriate for the task that needs to be completed.

The National Institute for Occupational Safety and Health (NIOSH) states the need for appropriate glove protection around infectious materials. Further, the CDC reports that the selection of glove type for non-surgical use is based on a number of factors, including the task that is to be performed, anticipated contact with chemicals and chemotherapeutic agents, latex sensitivity, sizing, and facility policies for creating a latex-free environment.

The CDC also advises that a single pair of gloves generally provides adequate barrier protection for contact with blood and bodily fluids during non-surgical patient care. However, there is considerable variability among gloves. Both the quality of the manufacturing process and the type of material influence their barrier effectiveness. While there is little difference in the barrier properties of unused intact gloves, studies have shown repeatedly that vinyl gloves have higher failure rates than latex or nitrile gloves when tested under simulated and actual clinical conditions. For this reason, either latex or nitrile gloves are preferable for clinical procedures that require manual dexterity and/or will involve more than brief patient contact.

In summary, the physical requirements of tasks performed in the hospital setting range from low to high according to the risk of exposure to infectious materials. Separating tasks into barrier durability needs can assist the decision making process on the level of barrier protection needed.

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A Comparison of Synthetic Gloving Materials

Glove Material	Acceptable Barrier Level	Elasticity	Puncture Resistance	Fit and Comfort	Chemical Resistance	Economy	Low Allergen Content	Environmental Impact	Recommended Use
 Nitrile	Low Medium High	Very Good	Excellent	Very Good	Excellent	Good	Very Good	Varied Does not decompose in landfills. Incineration produces mostly water and carbon dioxide.	Nitrile gloves are recommended as an alternative to NRL products, for conditions where high strength and chemical protection are required, and for individuals with latex allergies.
 Neoprene	Low Medium High	Excellent	Good	Excellent	Excellent	Good	Excellent	Varied Does not decompose in landfills. Incineration produces significant amounts of hydrochloric acid.	Neoprene gloves are recommended as an alternative to NRL because of their similarity to latex with respect to fit, feel, comfort, and barrier protection. They are also recommended for individuals with latex allergies.
 Vinyl	Low	Fair	Fair	Fair	Fair	Fair	Very Good	Poor Does not decompose in landfills. The plasticizers in vinyl are hormone mimics: if they leach out into the environment, they can have significant impact on wildlife. Incineration produces significant amounts of hydrochloric acid.	Vinyl use should be limited to short procedures in low-risk situations that do not involve exposure to blood, body fluids, or chemicals.

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Canadian General Standards Board delivers independent certification programs where products are assessed against performance and quality standards. Certification programs are a voluntary and impartial mechanism by which suppliers can demonstrate the value of their goods.

Company name and qualified products are displayed on the program list on-line at: pwgsc.gc.ca/cgsb/prgsrv/certprg/program/q020_ind-e.html

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- Ackerman, Frank and Massey, Rachel; Global Development and Environment Institute, Tufts University. The Economics of Phasing Out PVC. http://www.ase.tufts.edu/gdae/Pubs/rp/Economics_of_PVC_revised.pdf
- CDC. Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings 2007. <http://www.cdc.gov/ncidod/dhqp/pdf/guidelines/Isolation2007.pdf>

For additional information and samples, call Ansell Canada at 1 800 363-8340

® & TM are trademarks owned by Ansell Healthcare Products LLC or one of its affiliates. ©2009 All Rights Reserved.

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® & TM sont des marques déposées de Ansell Healthcare Products LLC ou une société affiliée. ©2009 Tous droits réservés.