EZMEND⁺™

Knee & Hip Rehabilitation Device

FASTER RECOVERY
FROM KNEE + HIP SURGERY

www.ez-mend.com

EZMEND KNEE AND HIP REHABILITATION DEVICE PATIENTS' GUIDE

Following knee and hip surgery, weakening of the muscles around these joints is a common occurrence. If you want to minimize stiffness of muscles, and increase mobility and muscle tone, get an **EZMEND Knee and Hip Rehabilitation Device**. This device may assist after knee and hip surgery, sports injuries and for nursing home patients.

EZMEND is engineered to improve circulation, increase your range of motion, strengthen muscles, diminish post-operative pain and reduce swelling.

Orthopedic surgeons and physiotherapists recommend **EZMEND** to help restore mobility and strength, to promote a faster and simpler recovery.

EZMEND can help as you begin to exercise your joints after surgery, as it is a non-weight bearing exercise device. It is important that you begin to exercise your joints following surgery, under your surgeon's supervision, to help restore range of motion and strength.

INSTRUCTIONS FOR KNEE & HIP RECOVERY:

These are guideline instructions that vary from patient to patient. It will be necessary to consult your doctor or physiotherapist about what degree of measurement and flexibility they would like to achieve. They will also be able to give you directions regarding the frequency of your exercise.

- 1) Sit lengthwise on your bed and position one end of the **EZMEND** under your thigh, parallel to your leg.

- 2) Wrap the soft strip around the instep of your foot.

- 3) Place your heel into the **EZMEND** cup and lay back into a horizontal position.

- 4) Grasp one end of the soft strip in each hand and slowly pull your leg along the **EZMEND** towards your chest. Hold your knee in this position for 5 – 10 seconds, and then straighten. Rest when fatigued.

- 5) Repeat this procedure following the directions and supervision of your doctor or physiotherapist.

INSTRUCTIONS FOR ADDITIONAL HIP EXERCISES:

These are guideline instructions that vary from patient to patient. It will be necessary to consult your doctor or physiotherapist about what degree of measurement and flexibility they would like to achieve. They will also be able to give you directions regarding the frequency of your exercise.

- 1) Sit lengthwise on your bed and position the **EZMEND** perpendicular below your leg on the bed.

- 2) Wrap the soft strip around the instep of your foot.

- 3) Place your heel in the **EZMEND** cup and lay back into a horizontal position.

- 4) Grasp one end of the soft strip in each hand and slowly pull your leg left to right along the **EZMEND**. Rest when fatigued.

- 5) Repeat this procedure following the directions and supervision of your doctor or physiotherapist.

After performing these exercises for a period of time, your doctor or physiotherapist may recommend you perform them while sitting on the side of your bed or in a chair with the **EZMEND** on the floor. Be sure not to modify your exercise routine without consulting your doctor.

Don't be alarmed if, initially, your exercises are uncomfortable. Continue with the **EZMEND** exercises under your doctor or physiotherapist's supervision and recovery time will decrease while your overall quality of life improves.

It is important that your doctor or physiotherapist monitors your progress with the **EZMEND** throughout the entire recovery period.

INSTRUCTIONS - CARE PROCEDURES

The **EZMEND** device may be cleaned using a household or commercial disinfectant. The soft strip is machine washable.

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© 2012 2330207 Ontario Inc.