



Medglove Powdered  
MG150 Series



Characteristics and Overview:

Medglove Lightly Powdered Latex Examination Gloves are manufactured to the highest level of quality standards, which comply with and exceed the requirements of the ASTM, American Society for Testing and Materials, and the FDA, Food and Drug Administration.

These premium quality gloves deliver the unique advantages of comfort, dexterity and barrier protection. Features Smooth silky finish provides superb tactile-sensitivity. Ergonomically designed to fit extremely well for maximum comfort. Powdered with USP grade edible cornstarch for easy donning. Beaded cuff guards against tearing.

**Sizes**

MG150XS	MG150S	MG150M	MG150L	MG150XL
X-Small	Small	Medium	Large	X-Large
100 ea/bx	100 ea/bx	100 ea/bx	100 ea/bx	100 ea/bx

Glove Specifications

	ASTM	Medgluv
Length (mm)	Min 220	240
Palm width (mm)	95 ± 10	95 ± 10
Cuff Thickness (mm)	N/A	0.10 (3.9 mil)
Palm Thickness (mm)	min. 0.05	0.12 (4.7 mil)
Finger Thickness (mm)	min. 0.05	0.14 (5.5 mil)
Tensile Strength (MPa)	min. 14	29
Ultimate Elongation (%)	min 500	750

Glove Attributes

Material	Latex
Sterility	Non-Sterile
Color	Beige
Sizes	X-Small, Small, Medium, Large, X-Large
Glove Interior	Lightly Powdered
Glove Exterior	Smooth
Packaging	100 gloves per box, 10 boxes per case

Quality Standards

Exceeds current ASTM 3578 standard for critical defects (AQL 2.5). AQL for critical defects is 1.5.

Manufactured in accordance with Quality System ISO 9001

Caution: *No sensitizing threshold to natural rubber latex has been established. Any individual with hypersensitivity to natural rubber latex should consult a physician before using any latex product..*

