

Sterigluv Sterile MGSE100 Series



Characteristics and Overview:

Sterigluv sterile latex singles are now the ultimate choice use of medical professionals in high demand environments like Obstetricians and Gynecologists, intensive care unit ICC/CCU, traumas, burn centers and wound care that requires a high quality low protein latex sterile glove that delivers maximum barrier protection. Small package is easy to open and less storage.

Sizes

MGSE100S	MGSE100M	MGSE100L	MGSE100XL
Small	Medium	Large	X-Large
100 ea/bx	100 ea/bx	100 ea/bx	100 ea/bx



Glove Specifications

	ASTM	Sterigluv Singles
Length (mm)	Min 220	240
Palm width (mm)	95 ± 10	95 ± 10
Cuff Thickness (mm)	N/A	0.12 (5.0 mil)
Palm Thickness (mm)	min. 0.09	0.14 (5.5 mil)
Finger Thickness (mm)	min. 0.05	0.15 (6.2 mil)
Tensile Strength (MPa)	min. 14	25
Ultimate Elongation (%)	min 400	750

Glove Attributes

Material	Latex
Sterility	Sterile
Color	Beige
Sizes	Small, Medium, Large, X-Large
Glove Interior	Powder-Free, Chlorinated
Glove Exterior	Textured
Packaging	Singles - 100 gloves per box, 4 boxes per case

Quality Standards

Exceeds current ASTM D3578 standard for critical defects (AQL 2.5). AQL for critical defects is 1.5.
Manufactured in accordance with Quality System ISO 9001.

**** Sterilization Location: Malaysia**

**** CAUTION:** No sensitizing threshold to natural rubber latex has been established. Any individual with hypersensitivity to natural rubber latex should consult a physician before using any latex product.