

## TECHNICAL DATA SHEET

### Cohesive Bandage Roll Items # A1061-A to A6061-W

#### Intended Use

Cohesive bandage rolls used in non-invasive procedures such as the covering of a dressing or to support a body part as a result of an injury. Used in health care facilities.

Characteristic	Specification			
Description	Self-adherent bandage roll			
Material Composition	Latex	Flexible non-woven cloth (60- 70%): 84% Polypropylene / 12.5% Spandex / 3.5% Natural Latex		
		Adhesive Latex (30 - 40%): 91% Natural Latex / 5% White Oil / 2% Antioxidants		
		Item Code		
		A1061-B, A2061-B, A4061-B, A1061-T, A2061-T, A3061-T, A4061-T, A6061-T, A2061-R, A4061-R, A1061-A, A2061-A, A3061-A, A4061-A, A6061-A, A2061-PPA, A2061-PPA-Purina, A4061-PPA, A2061-MIX		
	Latex-Free	Flexible non-woven cloth (75-90%): 69-84% Polypropylene / 6% Spandex / 10-25% Synthetic Glue		
		Item Code		
		A1061-TLF, A2061-ALF, A2061-TLF, A2061-TLFS, A3061-ALF, A3061-TLF, A4061-ALF, A4061-TLF, A4061-TLFB, A4061-TLFS, A6061-TLF, A6061-TLFB, A6061-TLFS		
Color	Item Code		Color	
	A1061-B, A2061-B, A4061-B		Blue	
	A1061-T, A1061-TLF, A2061-T, A2061-TLF, A2061-TLFS, A3061-T, A3061-TLF, A4061-T, A4061-TLF, A4061-TLFB, A4061-TLFS, A6061-T, A6061-TLF, A6061-TLFB, A6061-TLFS		Tan or Skin	
	A2061-R, A4061-R		Red	
	A2061-ALF, A3061-ALF, A4061-ALF		Assorted colors	
	A1061-A, A2061-A, A3061-A, A4061-A, A6061-A		6 colors	
	A2061-PPA, A2061-PPA-Purina, A4061-PPA		Paw printed assorted colors	
	A2061-MIX		Assorted designs: 8 Bones, 7 Hearts, 7 Stars, 7 Blue camouflage, 7 Purple camouflage	
Dimension	Size	Item Code	Width (cm)	Length (m)
	1" x 5 yds	A1061-A, A1061-B, A1061-T, A1061-TLF	2,5	4,5
	2" x 5 yds	A2061-A, A2061-ALF, A2061-B, A2061-MIX, A2061-PPA, A2061-PPA-Purina, A2061-R, A2061-T, A2061-TLF, A2061-TLFS	5,0	
	3" x 5 yds	A3061-A, A3061-ALF, A3061-T, A3061-TLF	7,5	
	4" x 5 yds	A4061-A, A4061-ALF, A4061-B, A4061-PPA, A4061-R, A4061-T, A4061-TLF, A4061-TLFB, A4061-TLFS	10,0	
	6" x 5 yds	A6061-A, A6061-T, A6061-TLF, A6061-TLFB, A6061-TLFS	15,0	

## TECHNICAL DATA SHEET

### Cohesive Bandage Roll Items # A1061-A to A6061-W

Characteristic	Specification			
Packaging	Item Code	Rolls/Box	Boxes/Case	Rolls/Case
	A1061-A, A1061-B, A1061-T, A1061-TLF	72	6	432
	A2061-A, A2061-ALF, A2061-B, A2061-MIX, A2061-PPA, A2061-PPA-Purina, A2061-R, A2061-T, A2061-TLF	36		216
	A3061-A, A3061-ALF, A3061-T, A3061-TLF	24		144
	A4061-A, A4061-ALF, A4061-B, A4061-PPA, A4061-R, A4061-T, A4061-TLF	18		108
	A6061-A, A6061-T, A6061-TLF	12		72
	Item Code	Rolls/Case		
	A4061-TLFB	100 (NLT 98)		
	A6061-TLFB	72 (NLT 71)		
	Item Code	Roll/Pouch	Pouches/Case	Rolls/Case
	A2061-TLFS	1	36	36
	A4061-TLFS		18	18
	A6061-TLFS		12	12
	A4061-TLFB, A6061-TLFB	Rolls are placed in a corrugated shipping box and separated from one another with cardboard inserts and layer separators		

Item Code	Description	Specific Attributes
A1061-A, A1061-B, A1061-T, A1061-TLF	AMD-Ritmed Cohesive Bandages, 1" x 5 yds (2,5 cm x 4,5 m), Individually Wrapped, 72 Rolls/Box, 6 Boxes/Case, 432 Rolls/Case	Non-Sterile
A2061-A, A2061-ALF, A2061-B, A2061-MIX, A2061-PPA, A2061-PPA-Purina, A2061-R, A2061-T, A2061-TLF*	AMD-Ritmed Cohesive Bandages, 2" x 5 yds (5 cm x 4,5 m), Individually Wrapped, 36 Rolls/Box, 6 Boxes/Case, 216 Rolls/Case	Non-Sterile
A2061-TLFS	AMD-Ritmed Cohesive Bandages, 2" x 5 yds (7,5 cm x 4,5 m), Individually Wrapped, 1 Roll/Pouch, 36 Pouches/Case, 36 Rolls/Case	Sterile
A3061-A, A3061-ALF, A3061-T, A3061-TLF	AMD-Ritmed Cohesive Bandages, 3" x 5 yds (7,5 cm x 4,5 m), Individually Wrapped, 24 Rolls/Box, 6 Boxes/Case, 144 Rolls/Case	Non-Sterile
A4061-A, A4061-ALF, A4061-B, A4061-PPA, A4061-R, A4061-T, A4061-TLF*, A4061-TLFB*	AMD-Ritmed Cohesive Bandages, 4" x 5 yds (10 cm x 4,5 m), Individually Wrapped, 18 Rolls/Box, 6 Boxes/Case, 108 Rolls/Case	Non-Sterile
A4061-TLFS	AMD-Ritmed Cohesive Bandages, 4" x 5 yds (10 cm x 4,5 m), Individually Wrapped, 1 Roll/Pouch, 18 Pouches/Case, 18 Rolls/Case	Sterile
A6061-A, A6061-T, A6061-TLF*, A6061-TLFB*	AMD-Ritmed Cohesive Bandages, 6" x 5 yds (15 cm x 4,5 m), Individually Wrapped, 12 Rolls/Box, 6 Boxes/Case, 72 Rolls/Case	Non-Sterile
A6061-TLFS	AMD-Ritmed Cohesive Bandages, 6" x 5 yds (15 cm x 4,5 m), Individually Wrapped, 1 Roll/Pouch, 12 Pouches/Case, 12 Rolls/Case	Sterile



**TECHNICAL DATA SHEET**  
**Cohesive Bandage Roll**  
**Items # A1061-A to A6061-W**

**Regulatory**

FDA Class 1, FDA Product Code FQM

Health Canada Class 1

**Storage Conditions**

Store under ambient temperature and dry conditions

**Country of Origin**

China except for item codes A2061-TLFS, A4061-TLFS and A6061-TLFS, Canada