

TECHNICAL DATA SHEET

Cotton Ball Items # 53843 to 83874

Intended Use

Cotton balls used to absorb biological fluids, control bleeding or protect organs and/or structure from abrasion, drying and contamination. For internal use (X-Ray element) or external use (no X-Ray element)

Characteristic	Specification	
Description	Balls, in the shape of a sphere or a cylinder, made of a non-woven fibrous material.	
Material	Ball & String	100% Cotton
	X-Ray Element	Barium Sulfate (BaSO ₄) & PVC (Base Material)
Color	Cotton Ball	White
	String	Blue
	X-Ray Element	Blue
String	Item Code	String
	53843, 53850, 53858, 83841, 83851, 83858	Non-Strung
	73868, 83842, 83846, 83850, 83854, 83867XR, 83868, 83870, 83872, 83874	Double Strung
X-Ray Element	Item Code	X-Ray Element
	53850, 83841, 83842, 83846, 83850, 83854, 83867XR, 83870, 83872, 83874	Yes
	53843, 53858, 73868, 83851, 83858, 83868	No
Size	Size (inch)	Item Code
	$\frac{1}{8}$	83867XR, 83868
	$\frac{1}{4}$	83842
	$\frac{3}{8}$	73868, 83870
	$\frac{1}{2}$	83846
	$\frac{3}{4}$	83872
	1	53850, 83850, 83851
	1 $\frac{1}{4}$	53843, 83841, 83874
	1 $\frac{1}{2}$	83854
	2	53858, 83858

TECHNICAL DATA SHEET

Cotton Ball Items # 53843 to 83874

Characteristic		Specification						
Packaging	Sterile products (Strung)	Item Code	Balls / Count Card	Count Cards / Glassine	Glassines / Pouch	Pouches / Case	Balls / Case	
		83842, 83846, 83850, 83854	5	1	1	20	100	
		83868, 83870, 83872, 83874				50	250	
		73868				200	1000	
		83867XR				240	1200	
	Sterile products (Non-strung)	Item Code	Balls / Glassine	Glassine / Pouch	Pouches / Case	Balls / Case		
		83841, 83851	5	1	20	100		
		Item Code	Balls / Pouch	Pouches / Case	Balls / Case			
		83858	5	20	100			
	Non-sterile products (Non-strung)	Item Code	Balls / Glassine	Glassines / Case	Balls / Case			
		53843, 53850	5	20	100			
		Item Code	Balls / Bag	Bag / Case	Balls / Case			
		53858	1000	1	1000			

Item Code	Description	Color	Specific Attributes
53843	AMD-Ritmed Cotton Balls, 1 ¼", 5 Balls/Glassine, 20 Glassines/Case, 100 Balls/Case	White	Non-Sterile
53850	AMD-Ritmed Cotton Balls, 1", 5 Balls/Glassine, 20 Glassines/Case, 100 Balls/Case	White	X-Ray Element, Non-Sterile
53858	AMD-Ritmed Cotton Balls, 2", 1000 Balls/Bag, 1 Bag/Case, 1000 Balls/Case	White	Non-Sterile
73868	AMD-Ritmed Cotton Balls, ¾", 5 Balls/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 200 Pouches/Case, 1000 Balls/Case	White	Double Strung, Sterile
83841	AMD-Ritmed Cotton Balls, 1 ¼", 5 Balls/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Balls/Case	White	X-Ray Element, Sterile
83842	AMD-Ritmed Cotton Balls, ¼", 5 Balls/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Balls/Case	White	Double Strung, X-Ray Element, Sterile
83846	AMD-Ritmed Cotton Balls, ½", 5 Balls/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Balls/Case	White	Double Strung, X-Ray Element, Sterile
83850	AMD-Ritmed Cotton Balls, 1", 5 Balls/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Balls/Case	White	Double Strung, X-Ray Element, Sterile
83851	AMD-Ritmed Cotton Balls, 1", 5 Balls/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Balls/Case	White	Sterile

TECHNICAL DATA SHEET

Cotton Ball Items # 53843 to 83874

Item Code	Description	Color	Specific Attributes
83854	AMD-Ritmed Cotton Balls, 1 ½", 5 Balls/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Balls/Case	White	Double Strung, X-Ray Element, Sterile
83858	AMD-Ritmed Cotton Balls, 2", 5 Balls/Pouch, 20 Pouches/Case, 100 Balls/Case	White	Sterile
83867XR	AMD-Ritmed Cotton Balls, 1/8", 5 Balls/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 240 Pouches/Case, 1200 Balls/Case	White	Double Strung, X-Ray Element, Sterile
83868	AMD-Ritmed Cotton Balls, 1/8", 5 Balls/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 50 Pouches/Case, 250 Balls/Case	White	Double Strung, Sterile
83870	AMD-Ritmed Cotton Balls, 3/8", 5 Balls/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 50 Pouches/Case, 250 Balls/Case	White	Double Strung, X-Ray Element, Sterile
83872	AMD-Ritmed Cotton Balls, 3/4", 5 Balls/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 50 Pouches/Case, 250 Balls/Case	White	Double Strung, X-Ray Element, Sterile
83874	AMD-Ritmed Cotton Balls, 1 ¼", 5 Balls/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 50 Pouches/Case, 250 Balls/Case	White	Double Strung, X-Ray Element, Sterile

Regulatory

FDA Class 1, FDA Product Code GDY (X-Ray Element) and FDA Product Code NAC (No X-Ray Element)

Health Canada Class 2

Storage Conditions

Store under ambient temperature and dry conditions

