

TECHNICAL DATA SHEET

Cylindrical Sponge Items # 53217 to 83250

Intended Use

Non-woven cotton rolls intended to be used to absorb biological fluids, with or without X-Ray element.

Characteristic	Specification					
Description	Soft fibrous non-woven material in the shape of a cylinder					
Material Composition	Sponge	100% non woven cotton, latex-free				
	X-Ray Element	Barium Sulfate (BaSO ₄) & PVC (Base Material)				
Color	Sponge	White				
	X-Ray Element	Blue				
X-Ray Element	Item codes 53258: No X-Ray Other item codes: Single X-Ray visible at least at one end					
Strung	Item Code			Strung	String Position	
	53217, 53217-S, 53219, 53221, 53221-S, 53223, 73212, 73218, 73220, 73222, 73224, 83250			None	Not Applicable	
	53243, 53243-S, 73224A, 83242, 83244			Double	End	
	53245, 53258, 83246, 83248				Center	
Dimension	Size	Item Code		Diameter (mm)	Length (mm)	
	3/8" x 1/4"	73212		9,5	8,0	
	3/8" x 3/8"	83250			14,0	
	3/8" x 3/4"	53217, 53217-S, 73218			25,5	
	3/8" x 1"	53219, 73220, 83242			36,5	
	3/8" x 1 1/2"	53221, 53221-S, 53243, 53243-S, 73222, 83244		12,5	36,5	
	1/2" x 1 1/2"	53223, 53245, 53258, 73224, 73224A, 83246			51,0	
	1/2" x 2"	83248				
Packaging (Sterile products)	Item Code	Sponges/Glassine	Glassine/ Pouch	Pouches/ Case	Sponges/ Case	
	73220, 73222, 73224	5	1	20	100	
	Item Code	Sponges/Hol der	Holder/ Glassine	Glassine/ Pouch	Pouches/ Case	Sponges/ Case
	73212, 73218	5	1	1	20	100
	83250				40	200
	Item Code	Sponges/Co unt Card	Count Card/ Glassine	Glassine/ Pouch	Pouches/ Case	Sponges/ Case
73224A, 83242, 83244, 83246, 83248	5	1	1	20	100	



TECHNICAL DATA SHEET

Cylindrical Sponge Items # 53217 to 83250

Characteristic	Specification				
Packaging (Non-Sterile products)	Item Code	Sponges/ Glassine	Glassines/Case		Sponges/Case
	53219, 53221, 53221-S, 53223	5	100		500
	Item Code	Sponges/Holder	Holder/Glassine	Glassines/Case	Sponges/Case
	53217, 53217-S	5	1	100	500
	Item Code	Sponges/ Count Card	Count Card/ Glassine	Glassines/ Case	Sponges/ Case
53243, 53243-S, 53245, 53258	5	1	100	500	

⁽¹⁾ Frays that are unacceptable are those that could come off

⁽²⁾ Quantity not less than 98% of Label Claim. NLT: Not Less Than

⁽³⁾ For information purposes only

Item Code	Description	Color	Specific Attributes
53217	AMD-Ritmed Cylindrical Sponge, 3/8" x 3/4", 5 Sponges/Holder, 1 Holder/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, X-Ray, Holder
53217-S	AMD-Ritmed Cylindrical Sponge, 3/8" x 3/4", 5 Sponges/Holder, 1 Holder/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Pre-Treated, Non-Sterile, X-Ray, Holder
53219	AMD-Ritmed Cylindrical Sponge, 3/8" x 1", 5 Sponges/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, X-Ray
53221	AMD-Ritmed Cylindrical Sponge, 3/8" x 1 1/2", 5 Sponges/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, X-Ray
53221-S	AMD-Ritmed Cylindrical Sponge, 3/8" x 1 1/2", 5 Sponges/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Pre-Treated, Non-Sterile, X-Ray
53223	AMD-Ritmed Cylindrical Sponge, 1/2" x 1 1/2", 5 Sponges/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, X-Ray
53243	AMD-Ritmed Cylindrical Sponge, 3/8" x 1 1/2", 5 Sponges/Count Card, 1 Count Card/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, X-Ray, End Strung
53243-S	AMD-Ritmed Cylindrical Sponge, 3/8" x 1 1/2", 5 Sponges/Count Card, 1 Count Card/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Pre-Treated, Non-Sterile, X-Ray, End Strung
53245	AMD-Ritmed Cylindrical Sponge, 1/2" x 1 1/2", 5 Sponges/Count Card, 1 Count Card/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, X-Ray, Center Strung
53258	AMD-Ritmed Cylindrical Sponge, 1/2" x 1 1/2", 5 Sponges/Count Card, 1 Count Card/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, Center Strung
73212	AMD-Ritmed Cylindrical Sponge, 3/8" x 1/4", 5 Sponges/Holder, 1 Holder/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, X-Ray, Holder



TECHNICAL DATA SHEET

Cylindrical Sponge Items # 53217 to 83250

Item Code	Description	Color	Specific Attributes
73218	AMD-Ritmed Cylindrical Sponge, 3/8" x 3/4", 5 Sponges/Holder, 1 Holder/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, X-Ray, Holder
73220	AMD-Ritmed Cylindrical Sponge, 3/8" x 1", 5 Sponges/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, X-Ray
73222	AMD-Ritmed Cylindrical Sponge, 3/8" x 1 1/2", 5 Sponges/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, X-Ray
73224	AMD-Ritmed Cylindrical Sponge, 1/2" x 1 1/2", 5 Sponges/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, X-Ray
73224A	AMD-Ritmed Cylindrical Sponge, 1/2" x 1 1/2", 5 Sponges/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, X-Ray, End Strung
83242	AMD-Ritmed Cylindrical Sponge, 3/8" x 1", 5 Sponges/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, X-Ray, End Strung
83244	AMD-Ritmed Cylindrical Sponge, 3/8" x 1 1/2", 5 Sponges/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, X-Ray, End Strung
83246	AMD-Ritmed Cylindrical Sponge, 1/2" x 1 1/2", 5 Sponges/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, X-Ray, Center Strung
83248	AMD-Ritmed Cylindrical Sponge, 1/2" x 2", 5 Sponges/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, X-Ray, End Strung
83250	AMD-Ritmed Cylindrical Sponge, 3/8" x 3/8", 5 Sponges/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 40 Pouches/Case, 200 Sponges/Case	White	Sterile, X-Ray, Holder

Regulatory

FDA Class 1, FDA Product Code GDY (X-Ray), FDA Product Code NAC (No X-Ray)
Health Canada Class 2, Canadian License Number 4354

Storage Conditions

Store under ambient temperature and dry conditions

Country of Origin

Canada