

TECHNICAL DATA SHEET

Procedure Earloop Face Mask, Level 2 Barrier Items # 2115, 2115K, 2116, 2117 and 2118

Intended Use

The procedure masks are intended to be worn by personnel during surgical or other medical procedures to protect both the patient and the operating personnel from transfer of microorganisms, body fluids and particulate material transfer.

Characteristic			Specification			
Description			Rectangular face masks with a shapeable nosepiece and two earloops present, one on each side, in order to hold mask in place.			
Material	Mask body Composition		Outer Material	Spunbond Polypropylene (SBPP)		
			Filter Layer	Meltblown Polypropylene (MBPP)		
			Inner Material	Spunbond Polypropylene (SBPP)		
	Nose-piece Composition		Aluminum / Plastic covered iron			
	Elastic Band Composition		Knitted polyester			
Color			Item Code	Color		
			2115	Blue (Pantone 9440C – 9442C)		
			2115K	White		
			2116	Pink (Pantone 9021C – 9286C)		
			2117	Yellow (Pantone 600C – 601C)		
			2118	Green (Pantone 331C – 333C)		
Dimension	Body Size	Length	Item 2115 K		Other Items	
			145 mm		175 mm	
		Width	95 mm		89 mm	
	Length of Earloop		160 mm		165 mm	
	Pleat Depth		14 mm			
	Length of Nose-piece		100 mm		105 mm	
Packaging			Item Code	Masks / Box	Boxes / Case	Masks / Case
			2115, 2115K, 2116, 2117, 2118	50	10	500
			Item Code	Masks / Plastic Bag	Plastic Bags / Case	Masks / Case
			2115-Bulk, 2117-Bulk	500	1	500
			2117-Fold	50	10	

TECHNICAL DATA SHEET

Procedure Earloop Face Mask, Level 2 Barrier Items # 2115, 2115K, 2116, 2117 and 2118

Visual Mask



Item Code	Description	Color	Specific Attributes
2115	Distech Procedure Earloop Face Masks, 50 Masks/Box, 10 Boxes/Case, 500 Masks/Case	Blue	Non-Sterile
2115-Bulk	Distech Procedure Earloop Face Masks, 500 Masks/Bag, 1 Bag/Case	Blue	Non-Sterile Bulk
2115K	Distech Procedure Earloop Face Masks, 50 Masks/Box, 10 Boxes/Case, 500 Masks/Case	White	Non-Sterile
2116	Distech Procedure Earloop Face Masks, 50 Masks/Box, 10 Boxes/Case, 500 Masks/Case	Pink	Non-Sterile
2117	Distech Procedure Earloop Face Masks, 50 Masks/Box, 10 Boxes/Case, 500 Masks/Case	Yellow	Non-Sterile
2117-Bulk	Distech Procedure Earloop Face Masks, 500 Masks/Bag, 1 Bag/Case	Yellow	Non-Sterile Bulk
2117-Fold	Distech Procedure Earloop Face Masks, 50 Masks/Plastic Bag, 10 Plastic Bags/Case, 500 Masks/Case	Yellow	Non-Sterile Mask folded in half
2118	Distech Procedure Earloop Face Masks, 50 Masks/Box, 10 Boxes/Case, 500 Masks/Case	Green	Non-Sterile

Regulatory

FDA Class 2, Product Code FXX, 510(K) number K051291

Health Canada Class 1

Storage Conditions

Store under ambient temperature and dry conditions. Avoid excessive heat (over 40°C or 104°F)