



## TECHNICAL DATA SHEET

### Surgical Tape Items # A5005 to A5230-B

#### Intended Use

Surgery tape used to fix medical devices during surgery.

Characteristic	Specification			
Description	White adhesive tape (plastic, paper or silk cloth) rolled on a plastic cylinder except item code A5018 which is rolled on a paper cylinder			
Material	Item Code	Fabric	Adhesive	
	A5005, A5005-B, A5010, A5010-18, A5010-B, A5018, A5020, A5020-B, A5034-18, A5010-15FT	Plastic	Polyethylene	
	A5105, A5105-B, A5110-B, A5110, A5120, A5120-B, A5130, A5130-B, A5110-15F	Paper	Polyester	
	A5205, A5205-B, A5210, A5210-B, A5220, A5220-B, A5230, A5230-B	Silk Cloth	Viscose & Polyester	
Color	Item Code		Color	
	A5005, A5005-B, A5010, A5010-18, A5010-B, A5018, A5020, A5020-B, A5034-18, A5010-15FT		Clear Matt	
	A5105, A5105-B, A5110-B, A5110, A5120, A5120-B, A5130, A5130-B, A5110-15F		Egg White	
	A5205, A5205-B, A5210, A5210-B, A5220, A5220-B, A5230, A5230-B		Egg White	
Dimension	Size	Item Code	Width (mm)	Length (m)
	(½" x 10 yds)	A5005, A5005-B, A5105, A5105-B, A5205, A5205-B	13	9,1
	(¾" x 18")	A5034-18	19	≥ 0,457
	(1" x 18")	A5018, A5010-18	25	
	(1" x 15')	A5010-15FT, A5110-15FT	25	4,6
	(1" x 10 yds)	A5010, A5010-B, A5110, A5110-B, A5210, A5210-B	25	9,1
	(2" x 10 yds)	A5020, A5020-B, A5120, A5120-B, A5220, A5220-B	50	9,1
(3" x 10 yds)	A5130, A5130-B, A5230, A5230-B	75	9,1	
Packaging <sup>(1)</sup>	Item Code	Rolls/Box	Boxes/Case	Rolls/Case
	A5130, A5130-B, A5230, A5230-B	4	12	28
	A5020, A5020-B, A5120, A5120-B, A5220, A5220-B	6		72
	A5010, A5010-B, A5110, A5110-B, A5210, A5210-B	12		144
	A5005, A5005-B, A5105, A5105-B, A5205, A5205-B	24	30	288
	A5010-18, A5034-18			720
	A5018	48	90	4320
	Item Code	Rolls/Case		
A5010-15FT, A5110-15FT	144			

<sup>(1)</sup> Items codes with a "B" at the end are sold per box



## TECHNICAL DATA SHEET

### Surgical Tape Items # A5005 to A5230-B

Item Code	Description	Color	Specific Attributes
A5005	AMD-Ritmed Surgical Tape; ½" X 10 yds (1,25 cm x 9 m), 24 Rolls/Box, 12 Boxes/Case, 288 Rolls/Case	Clear Matt	Plastic, Latex Free, Non-Sterile
A5005-B	AMD-Ritmed Surgical Tape; ½" X 10 yds (1,25 cm x 9 m), 24 Rolls/Box	Clear Matt	Plastic, Latex Free, Non-Sterile
A5010	AMD-Ritmed Surgical Tape, 1" X 10 yds (2,5 cm x 9 m), 12 Rolls/Box, 12 Boxes/Case, 144/Case	Clear Matt	Plastic, Latex Free, Non-Sterile
A5010-18	AMD-Ritmed Surgical Tape, 1" X 18" (2,5 cm x 0,456 m), 24 Rolls/Box, 30 Boxes/Case, 744 Rolls/Case	Clear Matt	Plastic, Latex Free, Non-Sterile
A5010-B	AMD-Ritmed Surgical Tape, 1" X 10 yds (2,5 cm x 9 m), 12 Rolls/Box	Clear Matt	Plastic, Latex Free, Non-Sterile
A5010-15FT	AMD-Ritmed Surgical Tape, 1" X 15' (2,5 cm x 0,456 m), 144 Rolls/Case	Clear Matt	Plastic, Latex Free, Non-Sterile
A5018	AMD-Ritmed Surgical Tape, 1" X 18' (2,5 cm x 4,57 m), 48 Rolls/Box, 90 Boxes/Case, 4320 Rolls/Case	Clear Matt	Plastic, Latex Free, Non-Sterile
A5020	AMD-Ritmed Surgical Tape, 2" X 10 yds (5 cm x 9 m), 6 Rolls/Box, 12 Boxes/Case, 72 Rolls/Case	Clear Matt	Plastic, Latex Free, Non-Sterile
A5020-B	AMD-Ritmed Surgical Tape, 2" X 10 yds (5 cm x 9 m), 6 Rolls/Box	Clear Matt	Plastic, Latex Free, Non-Sterile
A5034-18	AMD-Ritmed Surgical Tape, ¾" x 18" (1,9 cm x 0,456 m), 24 Rolls/Box, 30 Boxes/Case, 744 Rolls/Case	Clear Matt	Plastic, Latex Free, Non-Sterile
A5105	AMD-Ritmed Surgical Tape, ½" X 10 yds (1,25 cm x 9 m), 24 Rolls/Box, 12 Boxes/Case, 288 Rolls/Case	Egg white	Paper, Latex Free, Non-Sterile
A5105-B	AMD-Ritmed Surgical Tape, ½" X 10 yds (1,25 cm x 9 m), 24 Rolls/Box	Egg white	Paper, Latex Free, Non-Sterile
A5110	AMD-Ritmed Surgical Tape, 1" X 10 yds (2,5 cm x 9 m), 12 Rolls/Box, 12 Boxes/Case, 144 Rolls/Case	Egg white	Paper, Latex Free, Non-Sterile
A5110-B	AMD-Ritmed Surgical Tape, 1" X 10 yds (2,5 cm x 9 m), 12 Rolls/Box	Egg white	Paper, Latex Free, Non-Sterile
A5110-15FT	AMD-Ritmed Surgical Tape, 1" X 15' (2,5 cm x 0,456 m) 144 Rolls/Case	Egg white	Paper, Non-Sterile, Latex Free, Bulk
A5120	AMD-Ritmed Surgical Tape, 2" X 10 yds (5 cm x 9 m), 6 Rolls/Box, 12 Boxes/Case 72 Rolls/Case	Egg white	Paper, Latex Free, Non-Sterile



## TECHNICAL DATA SHEET

### Surgical Tape Items # A5005 to A5230-B

Item Code	Description	Color	Specific Attributes
A5120-B	AMD-Ritmed Surgical Tape, 2" X 10 yds (5 cm x 9 m), 6 Rolls/Box	Egg white	Paper, Latex Free, Non-Sterile
A5130	AMD-Ritmed Surgical Tape, 3" X 10 yds (5 cm x 9 m), 4 Rolls/Box, 12 Boxes/Case, 48 Rolls/Case	Egg white	Paper, Latex Free, Non-Sterile
A5130-B	AMD-Ritmed Surgical Tape, 3" X 10 yds (5 cm x 9 m), 4 Rolls/Box	Egg white	Paper, Latex Free, Non-Sterile
A5205	AMD-Ritmed Surgical Tape, ½" X 10 yds (1,25 cm x 9 m), 24 Rolls/Box, 12 Boxes/Case, 288 Rolls/Case	Egg white	Silk Cloth, Latex Free, Non-Sterile
A5205-B	AMD-Ritmed Surgical Tape, ½" X 10 yds (1,25 cm x 9 m), 24 Rolls/Box	Egg white	Silk Cloth, Latex Free, Non-Sterile
A5210	AMD-Ritmed Surgical Tape, 1"X10 yds " (2,5 cm x 9 m), 12 Rolls/Box, 12 Boxes/Case, 144 Rolls/Case	Egg white	Silk Cloth, Latex Free, Non-Sterile
A5210-B	AMD-Ritmed Surgical Tape, 1" X 10 yds (2,5 cm x 9 m), 12 Rolls/Box	Egg white	Silk Cloth, Latex Free, Non-Sterile
A5220	AMD-Ritmed Surgical Tape, 2" X 10 yds (5 cm x 9 m), 6 Rolls/Box, 12 Boxes/Case, 72 Rolls/Case	Egg white	Silk Cloth, Latex Free, Non-Sterile
A5220-B	AMD-Ritmed Surgical Tape, 2" X 10 yds (5 cm x 9 m), 6 Rolls/Box	Egg white	Silk Cloth, Latex Free, Non-Sterile
A5230	AMD-Ritmed Surgical Tape, 3" X 10 yds (5 cm x 9 m), 4 Rolls/Box, 12 Boxes/Case 48 Rolls/Case	Egg white	Silk Cloth, Latex Free, Non-Sterile
A5230-B	AMD-Ritmed Surgical Tape, 3" X 10 yds (5 cm x 9 m), 4 Rolls/Box	Egg white	Silk Cloth, Latex Free, Non-Sterile

#### Regulatory

FDA Class 1, FDA Product Code KGX

Health Canada Class 1

#### Storage Conditions

Store under ambient temperature and dry conditions

#### Country of origin

China