

ULTIMATE PEACE OF MIND FOR ANY SURGICAL NEED WITH AN ADVANCED ALLERGY PROTECTION SOLUTION

A chemical accelerator-free glove that earns the GAMMEX® name

The precedent-setting **GAMMEX® Non-Latex Sensitive** surgical glove has the natural fit and comfort you love in GAMMEX® gloves, so you can focus exclusively on your performance and not your skin sensitivities.

Key features include:

- A major innovation in non-latex, chemical accelerator-free glove formulation, **SENSOPRENE®** allows for superior sensitivity while maintaining glove strength
- The chemical accelerator-free formulation ensures comprehensive protection for healthcare workers with chemical allergies to accelerators
- GAMMEX® design and manufacturing guarantee the desirable fit, donning properties, and grip you expect
- The soft cream color you want in a surgical glove

The **GAMMEX® Non-Latex Sensitive** surgical glove goes even further in helping to protect healthcare workers and patients against additional causes of skin reactions by:

- * Utilizing a formulation that does not contain casein – or any other chemicals containing animal derivatives.

The **GAMMEX® Non-Latex Sensitive** surgical glove also delivers the proven protection and quality you expect in all GAMMEX® gloves:

- Meets or exceeds ASTM standards
- Delivers strong resistance to viral permeation, verified by extensive bacteriophage test results as per ASTM F1671-07
- Provides resistance to alcohol-based hand disinfectants

Advantages of SENSOPRENE® glove formulation

SENSOPRENE® is a new formulation, designed for the ultimate sensitivity and comfort. Free of natural rubber latex and chemical accelerators, it helps reduce the risk of both latex and chemical allergies while also maintaining glove strength.



Dimensional Specifications

Length
(typical average values, mm)

305

(Size 5.5–Size 9)

Thickness (Single)

Finger (typical average values, mm)	Palm (typical average values, mm)	Cuff
0.140	0.130	0.120

Glove Material and Formulation

Color	Material	Accelerator-Free	Free of Animal Derivatives
Cream	Neoprene; no natural rubber latex	Yes	Yes

Product Numbers - Beaded Cuff

SIZE (identified inside glove)	ASPEN (Ansell Standard Product Number)
5 1/2	20277255
6	20277260
6 1/2	20277265
7	20277270
7 1/2	20277275
8	20277280
8 1/2	20277285
9	20277290

50 pairs per dispenser box; 4 boxes per case

