

# Freezpoint™

Colder than the Can!

## MULTI-USE DISPOSABLE CRYOSURGICAL DEVICE



**A safer, easier, and colder device**, that delivers an improved, low cost per procedure, cryosurgical option. The Freezpoint™ device is 27% - 62% colder than common multi-use disposable systems. The slim tip design aids in clear visualization of the lesion for precise delivery of  $-89^{\circ}\text{C}$ , non-flammable, Nitrous Oxide therapy. The outcome is a very safe and clinically effective treatment with reduced tissue damage to surrounding healthy skin. Freezpoint is conveniently available in 8 gram (CRY8) and 16 gram (CRY16) devices.

### Indications Include:

- Skin Tags
- Verruca Vulgaris
- Verruca Plantaris
- Molluscum Contagiosum
- Lentigo
- Actinic Keratosis
- Verruca Plana
- Seborrheic Keratosis
- Genital Lesions

**Skin Tags:** Treatment time 5-10 seconds



BEFORE TREATMENT



POST TREATMENT



DAY 7



DAY 14