

Cutisorb® Ultra




**Super absorbent
wound dressing**

Cutisorb® Ultra – Superior exudate management that wo

Cutisorb® Ultra combines innovative materials and state of the art manufacturing processes to provide superior absorption that works under pressure in moderate to highly exuding wounds. The incorporated distribution layer insures exudate is channeled throughout the dressing where the super absorbent particles quickly absorb and irreversibly bind moisture. The soft wound contact layer and superior absorption capacity insure gentle and less frequent dressing changes for a better quality of life during treatment.

Innovative 5-layer approach to exudate management




Superior absorption and fluid management

Cutisorb® Ultra absorbs more than four times the exudate of other wound dressings and permanently locks it away from the wound and surrounding skin even under compression or pressure.

High capacity means less frequent dressing changes providing more undisturbed healing time and saving time, money and resources.

* Data on file






Secure fluid management

Cutisorb® Ultra channels wound exudate to the dressing core where it is irreversibly bound to the super absorbers. The protective outer layer works with the super absorbers to protect the wound and prevent strikethrough and leakage even under pressure.



- 1 Cutisorb® Ultra
- 2 Standard absorption dressing


Cutisorb® Ultra: See the difference, experience the benefits



Indications

Cutisorb® Ultra is indicated for the management of moderate to highly exuding wound including

- pressure ulcers
- diabetic foot ulcers
- venous leg ulcers
- infected wounds under medical supervision



Cutisorb® Ultra:

Superior exudate management that works under pressure

- highly absorbent
 - ▶ fewer dressing changes
- maintains optimal wound moisture levels
 - ▶ faster wound healing
- super absorber retains fluid under compression
 - ▶ protects against maceration
- low profile
 - ▶ minimizes pressure points under compression or over bony prominences
- soft and conformable
 - ▶ for maximum patient comfort





Wound Phase

-  Necrotic
-  Infected
-  Sloughy
-  Granulating
-  Epithelializing

Wound Depth

-  Shallow
-  Deep
-  Shallow + deep

Exudate

-  Dry to low
-  Low to medium
-  Medium to high
-  High to excessive

Cutisorb® Ultra

Ref. No.	Size (cm/inches)	Dressings/ Box	HCPCS
72639-00	10 x 10 cm 4" x 4"	20	A6251
72639-01	10 x 20 cm 4" x 8"	20	A6252
72639-02	20 x 20 cm 8" x 8"	20	A6252
72639-03	20 x 30 cm 8" x 12"	20	A6253

Suggested HCPCS Code

HCPCS Code listed is only suggested code; other codes may be more applicable depending on product use. Codes are subject to change without notice.

Please refer to the product label and /or package insert for full instructions on the safe use of these products.



BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com
To order toll-free: BSN medical (+1) 800 552 1157

