

BASX™ BLADELESS TROCAR

Striking an affordable balance between clinical requirements and economic demands, BASX™ trocars offer only the essential features needed for general procedures and a limited number of product codes for efficient inventory management.

Durable Universal Seal

- No reducer caps needed
- Accommodates a wide range of instrument sizes (4.7mm–12.9mm)



Integrated Stability Thread Design

Recessed Stopcock Valve

- Mitigates inadvertent opening/closing of valve during surgery



Low Profile Housing



**Ethicon
Endo-Surgery**

PART OF THE JOHNSON & JOHNSON FAMILY OF COMPANIES

PRODUCT FACT SHEET

BASX™ BLADELESS TROCAR

BASX™ Bladeless Trocar

CODE	DESCRIPTION	LENGTH
TB5ST	Stability sleeve, 5mm diameter	75mm
TB5LT	Stability sleeve, 5mm diameter	100mm
TB11LT	Stability sleeve, 11mm diameter	100mm
TB12ST	Stability sleeve, 12mm diameter	75mm
TB12LT	Stability sleeve, 12mm diameter	100mm

BASX™ Sleeves

CODE	DESCRIPTION	LENGTH
CTB5ST	Stability sleeve, 5mm diameter	75mm
CTB5LT	Stability sleeve, 5mm diameter	100mm
CTB11LT	Stability sleeve, 11mm diameter	100mm
CTB12LT	Stability sleeve, 12mm diameter	100mm

BASX™ trocars are available only with a nonoptical bladeless obturator.



**Ethicon
Endo-Surgery**

PART OF THE  FAMILY OF COMPANIES