


LIGACLIP®

Designed for **versatility**  
and great **value**

LIGACLIP® Multi-Patient Use Single Clip Appliers



## Features

- Lateral and transverse grooves on clip surface to help provide on-structure clip security
- Proximal to distal clip closure designed to provide minimal gap
- Grooves on the inside jaw surfaces of applicator designed to increase in-jaw clip security
- Tumbled tips to help reduce trauma to vessels
- Angled jaws designed to improve visibility beyond structures
- A mechanical stop in the LX100 & LX200 Series designed to limit the maximum amount of force applied to tissue.

**ETHICON**  
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Shaping  
the future  
of surgery

## Product Overview

- LIGACLIP® Multi-Patient Use Single Clip Applicators provide a complete line of conventional, angled and right-angle products.
- The V-shaped clip is designed for ligation of tubular structures.
- The handles and clip cartridges are conveniently color-coded.
- Right-angle applicators include the LX130, LX230, LC330, LC430 models.

## Ordering Information

Code	Description	Tip Angle Specification	Clip Size	# Clips	Clip Code
LX107	7-1/2", Small clips, Blue handle	15° angle	Small	6 20	LT100 <span style="color: blue;">■</span>
LX110	10-1/2", Small clips, Blue handle	15° angle			LT102 <span style="color: blue;">■</span>
LX130	10-1/2", right angled, Small clips, Blue handle	90° angle			
LX205	5-3/4", Medium clips, Silver handle	15° angle	Medium	6 20	LT200 <span style="color: gray;">■</span>
LX207	7-1/2", Medium clips, Silver handle	15° angle			LT202 <span style="color: gray;">■</span>
LX210	10-1/2", Medium clips, Silver handle	15° angle			
LX220	10-1/2", angled jaw, Medium clips, Silver handle	45° angle			
LX230	10-1/2", right angled, Medium clips, Silver handle	90° angle			
LC307	7-1/2", Medium/Large clips, Green handle	15° angle	Med/Large	6	LT300 <span style="color: green;">■</span>
LC310	10-1/2", Medium/Large clips, Green handle	15° angle			
LC320	10-1/2", angled jaw, Medium/Large clips, Green handle	15° angle			
LC407	7-1/2", Large clips, Yellow handle	15° angle	Large	6	LT400 <span style="color: yellow;">■</span>
LC410	10-1/2", Large clips, Yellow handle	15° angle			
LC420	10-1/2", angled jaw, Large clips, Yellow handle	45° angle			
LC430	10-1/2", right angled, Large clips, Yellow handle	90° angle			

Instruments per Box: 1