

Tubigrip™ Arthro-pad

FOAM PADDED TUBULAR ELASTICATED SUPPORT BANDAGE

- Tubigrip Arthro-pad provides non-invasive mechanical support for the elbow or knee after injury or surgery

- Washable and reusable

- Easy to apply and remove

- Provides firm, comfortable support for the affected joint

- Retains body heat



HOW TUBIGRIP ARTHRO-PAD WORKS

Tubigrip Arthro-pads are constructed from polyurethane foam enclosed in a double layer of Tubigrip tubular bandage.

This double layer of Tubigrip provides firm support, encouraging use of the affected joint. In addition, the combination of Tubigrip and polyurethane foam retains body heat.

Arthro-pads are quick and easy to apply, making them suitable for outpatient and community use.

They are also washable, for repeated use. The neutral colour of Tubigrip Arthro-pad is cosmetically acceptable to most patients.

BENEFITS OF TUBIGRIP ARTHRO-PAD

- Provides firm, comfortable support for the affected joint
- Retains body heat
- Easily applied and removed
- Washable, for repeated use

MATERIAL PROPERTIES

- Polyurethane foam
- Elasticated support bandage consisting of:
 - Cotton
 - Elastodiene (contains natural rubber latex)
 - Polyamide

Washing instructions



AREAS OF USE

Tubigrip Arthro-pads are used to provide support for the elbow or knee after injury or surgery. They may also be of value to patients with arthritis of these joints.

TUBIGRIP ARTHRO-PAD ORDERING INFORMATION

Art. No	Size	Joint Circumference	Pieces Per Case	NHS Code
1585	Small	Up to 30cm	40	ESM011
1586	Medium	Up to 40cm	40	ESM013
1587	Large	Up to 50cm	40	ESM015
1588	X-Large	Up to 60cm	40	ESM017
1589	XX-Large	Up to 70cm	40	ESM017

Supplied individually packed



Mölnlycke Health Care Ltd, Arenson Centre, Arenson Way, Dunstable, Bedfordshire LU5 5UL, UK. Tel. +44 (0) 870 6060766, www.molnlycke.co.uk

The Mölnlycke Health Care name and logo, Tubigrip® Arthro-Pad are registered trademarks of Mölnlycke Health Care AB.
©2012 Mölnlycke Health Care. All rights reserved.

UK139550611

