

800-365-3047 • 715-294-2050
www.coreproducts.com

- Exceptional ankle stabilization and arch support
- Orthotic foot plate helps prevent inversion and eversion
- Provides stabilization without disrupting normal gait
- Neutral arch support promotes healing of injured ligaments



Sized by US shoe size

RIGHT MODEL #	LEFT MODEL #	SIZE	MEN'S	WOMEN'S
AKL-6372-R-BLK-1XS	AKL-6372-L-BLK-1XS	XS	-	4-6
AKL-6372-R-BLK-SML	AKL-6372-L-BLK-SML	SM	4-7	7-9
AKL-6372-R-BLK-MED	AKL-6372-L-BLK-MED	MED	8-10	10-12
AKL-6372-R-BLK-LRG	AKL-6372-L-BLK-LRG	LRG	11-13	13-15
AKL-6372-R-BLK-1XL	AKL-6372-L-BLK-1XL	XL	14+	16+



Critical Arch & Foot Support

The Swede-O Arch Lok® provides critical ankle and arch support. The modified orthotic footplate and ankle uprights help prevent inversion and eversion injuries without significant reduction of normal plantar and dorsiflexion so there is no negative impact to movement or performance.

This brace fits comfortably over the ankle and is secured with three adjustable straps. Once the brace is fit, easily apply and remove the brace using the single back strap. Use this product for patients with mild to moderate PTTD (Posterior Tibial Tendon Dysfunction) and ankle instability. Specify right or left.

Indications for Use

The Swede-O Arch Lok® provides exceptional ankle stabilization and arch support without disrupting norm gait. The orthotic foot plate helps prevent inversion and eversion of the foot.

Specifications

Product Selection:

#AKL-6372 Swede-O Arch Lok®

Product Weight: .5 lb

Dimensions: 10" x 13.75" x 4"

Case Weight: 4 lbs

Case Dimensions: 11" x 14" x 13"

Case Quantity: 6

Available From:

