

Electrotherapy

InTENSity™ Micro Combo

The InTENSity™ Micro Combo is a Microcurrent device that lets the patient switch from TENS therapy to Microcurrent therapy, allowing them treatment of a variety of musculoskeletal issues.



Product Features:

- Simple to operate, portable device for pain management
- Combines TENS and Microcurrent allowing patients to treat a wide variety of issues
- Comes complete with carrying case, two sets of lead wires, 9-volt battery, AC adaptor, manual and pack of four self-adhesive reusable electrodes
- Power Source: Standard 9-volt battery or AC adaptor
- Dimensions: 4.5" x 2.55" x 0.9"
- Weight: 0.28 lb. with battery
- One year warranty
- HCPCS Code: E0730* for TENS Unit
- HCPCS Code: A4595* for electrotherapy supplies

Product Specifications:

	TENS	Microcurrent
Wave Form	Monophasic Square	Monophasic Square
Pulse Rate	1-150 Hz	1-150 Hz
Pulse Width	50-300 μ s	120-200 μ s
Output Voltage	Constant Current	Constant Current
Current	0-105 mA	0-700 μ A
Modes	4	3
Timer	1-60 Constant	1-60 Constant

Item #	Description	Qty
DI9698	InTENSity™ Micro Combo	each

* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:

(800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com