


## Electrotherapy

# TENS 7000 To Go™ Back Pain Relief System

The TENS 7000 To Go™ Back Pain Relief System is a doctor recommended, over-the-counter (OTC) TENS device that provides prescription-strength pain relief associated with sore and aching muscles in the lower back. This easy-to-use device comes equipped with a conductive electrode support belt that transmits the current directly to the pain site. With five preset therapy modes, you can quickly customize your therapy and target your back pain. Comes complete with TENS 7000™ device, conductive electrode support belt, conductive spray, snap-on lead wire, a 9-volt battery and user manual. FDA cleared for over-the-counter sale, a great retail item.



### On-the-Spot Pain Relief



### Product Features:

- FDA cleared for over-the-counter sale (no prescription required)
- TENS technology stimulates body's own natural pain killers
- Simple-to-operate, portable device and support belt for targeted pain management
- Device has five pre-set therapy modes allowing you to customize your therapy and target your pain
- Comes complete with TENS 7000™ device, conductive electrode support belt, conductive spray, snap-on lead wire, a 9-volt battery and user manual

### DRUG-FREE RELIEF for muscular aches & pain



### Technical Specifications:

- Channel: Dual channels
- Modes: Burst mode, normal mode, modulation mode, SD1 mode, SD2 mode
- Output: 0 - 80 mA (500 ohm)
- Treatment Time: 5 - 60 minutes or continuous mode
- Pulse Rate: Adjustable from 2 to 150 Hz, 1hz/step
- Pulse Width: Adjustable from 50 to 300 microseconds, 10µs/step
- Dimensions: 4" x 6.25" x 1" with out belt
- Device Weight: 0.33 lb. with battery and belt
- Power: Supply: 9-volt battery
- Warranty: One year warranty from date of purchase

Item #	Description	Qty
DT6070	TENS 7000 To Go™ Back Pain Relief System	each
GB7000	TENS 7000 To Go™ - Replacement Conductive Back Brace	each
LS2177	Conductive Spray, 2 ounce	each

\*This device is FDA cleared for over-the-counter sale for pain relief associated with sore and aching muscles in the lower back.



**To Order:**

(800) 3-ROSCOE  
fax: (440) 572-4261  
roscoemedical.com